

A-5035


Sihari Raghavendra Rao¹, Chandrashekhar Honrao¹, Keri Sheehan¹, Leo Rodrigues¹, Chen Dong¹, Elizabeth O'Day^{1*}, and Dirk R. Kuypers²
¹Olaris, Inc., Framingham, MA 01702, U.S.A. ²Dept. of Nephrology and Renal Transplantation, University Hospitals, Leuven, Belgium *Corresponding author

WHAT IS THE PROBLEM?


Managing complications related to over-immunosuppression represents a growing challenge in post-transplant care, with **infections accounting for the second leading cause of death with functioning graft (DWFG)** in renal transplant recipients (RTRs) within the first year¹. At present, there are no clinically validated biomarkers to detect over-immunosuppression². **Polyomavirus-associated nephropathy (PVAN)** is a specific type of opportunistic infection indicative of over-immunosuppression that occurs in 5-10% of RTRs and causes up to 60% of graft failures³. Biopsy remains the gold standard for PVAN diagnosis, and only strongly increased blood viral copy number is associated with an increased risk for PVAN. Using urine, we leveraged metabolic profiling and machine learning to develop a novel metabolite-based signature to differentiate RTRs with a stable graft (no rejection or infection) from those developing biopsy-confirmed PVAN during follow-up. These results could lead to a non-invasive assay for the early detection of PVAN.

THE OLARIS CEREBRO PLATFORM COMBINES METABOLOMICS & MACHINE LEARNING

Altered metabolism has been linked to kidney function and kidney transplant tolerance. By measuring the complete set of metabolites in an individual (their metabolome), it is possible to identify biomarkers that correlate with disease status, prognosis, and therapeutic response. We detect and quantify **metabolites** from hundreds of patient biofluid samples (blood and/or urine). Using **NMR**- and **MS**- based metabolomics and **machine learning**, we identify a biomarker-based score that can monitor patient response to treatment with a high degree of accuracy.




STUDY DESIGN



Urine metabolites were extracted and analyzed via 1D ¹H and 2D ¹H-¹³C Heteronuclear Spectrum Quantum Coherence (HSQC) NMR spectroscopy of samples from **70 RTRs collected 3 months post-transplant, including N=38 with subsequent biopsy-proven PVAN (labeled "Over") and N=32 with stable graft function and no histological signs of rejection or indications of over-immunosuppression for 2 years (labeled "Stable")**. Spectra were processed using in-house processing and normalization tools and metabolite resonance peaks were assigned to metabolites based on chemical shift mapping to a library of known metabolites. Differential metabolite resonances between Over and Stable samples were then used to generate cross-validated machine learning models, including orthogonal partial least square discriminate analysis (OPLS-DA) and random forest (RF). The model with the highest stability and overall performance was selected as the champion model. Based on the champion model, we then developed a scoring method (**myOLARIS-KTdx™ Score**) with high accuracy to differentiate Stable and Over RTRs.

DIFFERENTIAL METABOLITE RESONANCES IN OVER VS. STABLE RTRs




Annotation	1H ppm	13C ppm	Fold change (Stable/Over)	P-value	FDR
Trigonelline	9.117	148.323	1.74	0.001	0.052
Trigonelline	8.841	147.243	1.578	0.001	0.052
Trigonelline	8.08	130.247	1.567	0.006	0.087
4-Ethylbenzoic acid	7.817	131.399	2.405	0.007	0.087
Unknown_7_531_129.748	7.53	129.739	2.268	0.001	0.056
Unknown_7_532_134.804	7.531	134.794	2.566	0.003	0.065
3-Hydroxymandelic acid,L-Phenylalanine,4-Methoxyphenylacetic acid	7.292	132.713	1.758	0.007	0.092
Unknown_7_346_121.584	7.347	121.584	1.73	0.006	0.087
Unknown_7_085_121.7	7.084	121.706	2.047	0.002	0.058
Salicylic acid	6.974	120.638	1.729	0.006	0.087
Unknown_7_274_116.531	7.273	116.5	1.967	0.002	0.056
Homoveratric acid	6.914	115.412	1.758	0.005	0.084
Quinic acid	1.994	40.118	1.602	0.002	0.058
Hippuric acid	7.818	129.759	1.701	0.006	0.087
Trigonelline	4.441	50.943	1.627	0.001	0.052
Unknown_3_813_62.593	3.81	62.602	1.509	0.024	0.19
Unknown_2_712_47.752	2.712	47.748	0.598	0.038	0.261
Unknown_2_566_47.724	2.566	47.755	0.569	0.022	0.185
Unknown_3_787_73.579	3.785	73.587	0.655	0.044	0.293
4-Aminohippuric acid, Phenylglylamic acid	7.624	131.127	1.646	0.008	0.095
Hippuric acid	7.618	134.806	2.494	0.001	0.052
Hippuric acid	7.531	131.362	2.33	0.003	0.056
Hippuric acid	3.969	46.485	2.04	0.003	0.056
Unknown_3_876_35.932	3.874	35.924	2.107	0.017	0.163

Using a Kruskal-Wallis (KW) non-parametric one-way analysis of variance, we identified **24 differential metabolite resonances** with fold-change >1.5 and p-value <0.05 between Over and Stable patients. The false discovery rate (FDR) was calculated to correct for multiple hypothesis testing, with many metabolite resonances with an FDR <0.1. Using chemical shift mapping, many differential resonances were mapped to the same metabolite. Resonances requiring experimental validation to distinguish matching-ambiguities are listed as "Unknown" and tracked by chemical shift.

GENERATING THE CHAMPION MODEL


Known Metabolic Resonances	Unknown Metabolic Resonances	Clinical Features
3-Hydroxymandelic acid,L-Phenylalanine, 4-Methoxyphenylacetic acid, 7.292_132.713, 4-Aminohippuric acid, Pteroylglutamic acid, 7.624_131.127, 4-Ethylbenzoic acid, 7.817_131.399, Hippuric acid, 7.531_131.362, Hippuric acid, 7.618_134.806, Hippuric acid, 7.818_129.759, Hippuric acid, 7.531_131.362, Hippuric acid, 3.969_46.485, Salicylic acid, 6.974_120.638, Trigonelline, 4.441_50.943, Trigonelline, 8.08_130.247, Trigonelline, 9.117_148.323	Unknown_3_787_73.579, Unknown_2_566_47.724, Unknown_2_712_47.752, Unknown_3_813_62.593, Unknown_3_876_35.932, Unknown_7_085_121.7, Unknown_7_274_116.531, Unknown_7_346_121.584, Unknown_7_531_129.748, Unknown_7_532_134.804	gender, median_dose_mmf, leukopenia_episodes, mean_dose_mmf, num_HLA_A_locus, iqr75_BW, num_HLA_B_locus, iqr25_BW, num_HLA_DR_locus, median_BW, acute_rejection_episodes, mean_BW, iqr75_dose_steroids, starting_dose_tac, iqr25_dose_steroids, iqr75_dose_tac, median_dose_steroids, iqr25_dose_tac, mean_dose_steroids, median_dose_tac, mean_dose_tac, iqr75_level_tac, iqr25_level_tac, median_level_tac, mean_level_tac

Pre-defined feature sets were generated including **"Known Metabolic Resonances"** (differential metabolite resonances assigned to a unique metabolite), **"Unknown Metabolic Resonances"** (differential metabolite resonances requiring experimental clarification), and **"Clinical Features"** (available sample-donor information such as gender, leukopenia, HLA matching, drug dose, trough levels, etc.).




Data were split 70/30 into training and test data, 10x cross-validated, and assessed based on accuracy, precision, sensitivity, and specificity. **A total of 16 ML models (8 OPLS-DA and 8 RF) were generated using different combinations of "Known", "Unknown" and "Clinical" features** with and without Boruta feature selection. Boruta selects features based on the RF importance of the unchanged vs. shuffled feature and can be used as a feature reduction method⁴. The RF Model with "Known" metabolites and optimized with Boruta feature selection was selected as the champion model (red box) based on optimal performance in training/test cohorts across all metrics.

PERFORMANCE OF myOLARIS SCORE FOR DETECTION OF OVER-IMMUNOSUPPRESSED & STABLE RTRs




The champion model was converted into a proprietary **myOLARIS Score** that ranges from 0 to 1, such that Over patients (orange) have higher scores than Stable patients (blue). Applying the myOLARIS Score to the full dataset demonstrated **excellent separation between Over and Stable RTRs with an overall accuracy of 83%**. This suggests that a urine metabolite signature can identify patients who are at risk for developing PVAN from those that are at low risk for either a rejection or infection event, i.e., the stable patients. The high positive predictive value (PPV) of 81% and high negative predictive value (NPV) of 86% offer the exciting possibility that this assay could be used for both **rule-in early diagnosis of PVAN and rule-out diagnosis of any event (rejection or infection), potentially eliminating the need for a biopsy**.



myOLARIS-KTdx Score
 84% AUC
 83% Accuracy
 89% Sensitivity
 75% Specificity
 81% PPV
 86% NPV

CONCLUSIONS AND NEXT STEPS

Determining the appropriate balance of immunosuppressant remains a clinical challenge for kidney transplant, as well as all solid organ transplant. Underimmunosuppression can lead to acute and/or chronic graft rejection, while overimmunosuppression can lead to complications such as infection, cardiovascular disease, post-transplant diabetes and malignancy. While numerous tools are available to detect risk of under-immunosuppression, there are few, if any, tools to detect and prevent complications of over-immunosuppression. In this study, we demonstrated that a metabolite-based urine signature was able to classify over-immunosuppressed RTRs with PVAN from those with a stable graft with high accuracy. Validation studies in larger cohorts with longitudinal samples are underway, which could lead to a **powerful new biomarker for post-transplant monitoring, empowering early diagnosis of PVAN and providing physicians and patients confidence when appropriate immunosuppression has been achieved for stable graft function**.



REFERENCES ~ ACKNOWLEDGMENTS

- Hariharan S, Israni AK, Danovitch G. Long-Term Survival after Kidney Transplantation. N Engl J Med. 2021 Aug 19;385(8):729-743. doi: 10.1056/NEJMra2014530.
- Lamarque C, Orio J, Collette S, et al. BK Polyomavirus and the Transplanted Kidney: Immunopathology and Therapeutic Approaches. Transplantation. 2016 Nov;100(11):2276-2287. DOI: 10.1097/tp.0000000000001333.
- Management of polyomavirus-associated nephropathy in renal transplant recipients. Dirk R J Kuypers. Review Nat Rev Nephrol. 2012 Apr 17;8(7):390-402. doi: 10.1038/nrneph.2012.64.
- Kursa M.B., Jankowski A., Rudnicki W.R. Boruta—A System for Feature Selection. Fundam. Inform. 2010;101:271-285.

Acknowledgements: We would like to acknowledge the Institute for Life Sciences at University of Massachusetts Amherst and Dr. Jasna Fejzo for assistance with processing urine samples using NMR.