


V. BIGAND¹, F. BRUNEL¹, C. LAMY¹ and F. TRANCHEPAIN¹ 1 Laboratoires VIVACY, 252 rue Douglas Engelbart, Archamps Technopôle, 74160 Archamps, France v.bigand@vivacy.fr

INTRODUCTION



- Natural polysaccharide, widely used in the biomedical field
- Features : Biocompatibility, biodegradability and biological properties




Hyaluronic acid

- Disaccharide unit
- In aesthetic medicine :
- \rightarrow Used as dermal filler available in prefilled syringes.
- \rightarrow Injected into the skin tissue in order to smooth and/or lift the skin

Chemical modification to obtain an hydrogel

Crosslinking reaction = Creation of covalent bonds between the polymer chains → Formation of a 3D network with improved mechanical properties and longer persistence under the skin

Crosslinking reaction between HA and BDDE (Butanediol diglycidyl ether):



IPN-like technology: Vivacy's patented crosslinking method, consisting in two interpenetrated networks Stylage[®] range of crosslinked HA: Hydromax, S, M, L, XL and XXL

	DERMAL FILLING ECOSYSTEM		
Products class	Need	Injection plans	Tiss featu
Mesotherapy	Skin Hydration & Rejuvenation	Epidermis Dermo-epidermic junction Papillar dermis	Stiff bu tissu
Fillers	Wrinkle & Skin depression filling	Papillar & Reticular dermis	Dense t
Volumizers	Volume restoration	Hypodermis & Subcutaneous tissue	Loose t
	class Mesotherapy Fillers	Products classNeedMesotherapySkin Hydration & RejuvenationFillersWrinkle & Skin depression filling	Products classNeedInjection plansMesotherapySkin Hydration & RejuvenationEpidermis Dermo-epidermic junction Papillar dermisFillersWrinkle & Skin depression fillingPapillar & Reticular dermisVolume restorationHypodermis &

AIM

In order to achieve their clinical performances, dermal filling products have to fulfil three main functions: Injectability, persistence on site and mechanical support.

The objective of this study is thus to investigate the injectability and viscoelastic properties of the Stylage[®] range of products and emphasize their relevance with regard to their targeted clinical effect.

METHOD


<u>Rheological measurements</u> were carried out using a DHR-2 controlled stress rheometer (TA instruments) with a coneplate geometry.

Elastic modulus (G') and **loss factor** (tan δ) \rightarrow frequency sweep 0.1 to 5 Hz & Strain 0.8% **Yield stress** (τ_c): stress value at the moduli cross-over point \rightarrow strain sweep from 0.1 to 1000% strain at 1 Hz

Injectability test: A needle was connected to the syringe containing the product and the injection profile was assessed using a traction bench (Mecmesin). A set of 3 compression speeds (13, 50 and 100 mm/minute) was successively applied to the piston rod and the resulting forces to extrude the product were measured.


11TH WORLD BIOMATERIALS CONGRESS

CLINICAL PERFORMANCES AND INJECTION PLANS




support the tissue

1- Characterization related to the injection procedure (example with Stylage[®] M)



Elastic Modulus (G')


- Characterize the rigidity/elasticity of the product



✓ Gradual increase of the rigidity within the range

Meso: Low rigidity for an homogeneous distribution **Fillers**: Rigidity increasing with the injection depth Volumizers: High rigidity required to lift up overlying tissues and to restore volumes

✓ Smooth and comfortable injection for both the patient and the practitioner ✓ Correlation between yield stress and rigidity: Both parameters participate to the *limitation of spreading and to the filling/projection capacity*





Fillers: Medium smoothness to reach a compromise between adaptability and filling performance


Volumizers: Lower smoothness to optimize the volumizing effect

CONCLUSIONS

Global rheological characteristics of products in accordance with the different clinical indications


Mesotherapy Injection plan: close to the skin surface

Product requirements: High **adaptability** to avoid visible irregularities



Stylage XI






✓ Identical rheological profile before and after injection → The product recovers 100% of its mechanical properties immediately after injection


Yield stress (T_c)

- Spreading capacity : can be related to viscosity features
- Characterize the resistance to flow
- Capacity to remain on the injection site



✓ **Increase of the spreading resistance** within the range

Meso: High spreading capacity to allow product nappage within the skin **Fillers**: Spreading decreases along with the density of surrounding tissues **Volumizers**: Low spreading for a high retention on site within loose subcutaneous tissues



∠ Bic

MBC2020

pocter."

Virtual Environments Communication