

EDTA IMPROVES THE CELL COUNT RESULTS IN THE LABORATORY DIAGNOSIS OF PD PERITONITIS

Fatima Abdelaal¹, Jodie Earp², Karen Simms¹, Lavanya Kamesh¹, Miruna David².

Peritoneal Dialysis Unit¹, Department of Clinical Microbiology², Queen Elizabeth Hospital Birmingham, University Hospitals Birmingham, Birmingham, UK.

fatima.abdelaal@uhb.nhs.uk.



INTRODUCTION

- Peritonitis remains a common complication of peritoneal dialysis (PD) and is associated with significant morbidity, catheter loss and transfer to haemodialysis.
- Patients with peritonitis usually present with abdominal pain and cloudy peritoneal effluent.
- A presumptive diagnosis is made if the peritoneal fluid yields a white cell count (WCC) $>100/\text{mm}^3$ and the percentage of neutrophils is $>50\%$. Diagnosis may be made in the presence of $>50\%$ polymorphonuclear (PMN) cells independent of the absolute WCC among APD patients. Diagnosis of peritonitis is confirmed by a positive dialysate culture [1].
- Laboratory analysis can provide a total as well as a differential cell count (DCC), gram stain and culture for diagnosis of peritonitis from PD fluid sent in a universal container.
- Frequently PD fluids received in universal containers are clotted and therefore unsuitable for a cell count.
- Time delays in processing samples may lead to inaccurate results, since white cells degrade with time.
- Guidelines already recommend submitting ascitic fluid samples in ethylenediaminetetraacetic acid (EDTA) in order to perform cell count for the diagnosis of spontaneous bacterial peritonitis [2].
- A previous study has shown that leukocytes stability in dialysis effluents is increased when samples are drawn in tubes containing EDTA [3].

OBJECTIVES

- To compare the PD fluid DCC collected in universal containers and EDTA tubes.

METHODS


- The study was carried out between June 2013 and January 2014.
- Twenty PD fluid samples were collected in duplicate, both in a sterile universal container and an EDTA tube.
- 0.1% toluidine blue stain was added to 2-3 drops of the PD fluid sample to differentiate between neutrophils and other white cells.
- The cell counts were performed simultaneously on both types of specimens to provide a true comparison.
- Time delay in processing the specimens was recorded to evaluate its impact on the results.
- Gram stain and culture were only carried out on the centrifuged sterile universal container samples.
- Training of laboratory staff was straightforward.



RESULTS

- Of the 20 samples collected, 50% of specimens in the universal containers had a WCC of $>100/\text{mm}^3$.
- 4 out of the 20 samples (20%) were clotted and the cell count could be performed only from the EDTA specimen.
- Three out of the 4 showed a WCC of $>100/\text{mm}^3$.
- When the remaining (unclotted) 16 samples were compared with each counterpart, there was an average 14% increase in the PMN count in the EDTA compared to the universal sample.
- The PMN counts by the two methods showed a good correlation with an R^2 value of 0.88 ($p < 0.001$).

Comparison of peritoneal fluid PMN cell counts done on plain tube and EDTA tube samples (n=16)



CONCLUSION

- Out of all PD fluid samples submitted in universal containers for cell count in cases of suspected peritonitis, a significant proportion clot.
- Using specimen containers with EDTA may improve the diagnosis of PD peritonitis.
- This practice therefore may have an impact on the management of culture negative peritonitis.

ACKNOWLEDGEMENT

We would like to acknowledge Dr Martin Gill for providing the statistical analysis for this study.

REFERENCES

- Li, PK et al. International Society for Peritoneal Dialysis Peritoneal dialysis-related infections recommendations: 2010 update. *Perit Dial Int* 2010; 30 (4): 393-423.
- Moore KP, Aithal GP. Guidelines on the management of ascites in cirrhosis. *Gut* 2006;55 Suppl 6:vi1-vi12.
- Antonsen S, Pedersen FB, Wang P and the Danish Study Group on Peritonitis in Dialysis (DASPID). Leukocytes in peritoneal dialysis effluents. *Perit Dial Int* 1991; 11 (1): 43-47.

