PROTEIN-CREATININE RATIO IN 24-HOUR URINE VERSUS SPOT URINE IN THE FOLLOW-UP OF PROTEINURIC CKD

A. Shabaka, M. Garbiras, V. Lopez de la Manzanara, J. Delgado, MC Cardenas, L. Martin, M. Calvo, F. Hadad, I Ubeda, JA. Herrero Nephrology Department. Hospital Universitario Clinico San Carlos, Madrid.

INTRODUCTION

Although the gold standard for the quantification of proteinuria is the 24-hour sample urine protein, due to its numerous limitations (inconvenience and errors in its collection), measurement of proteincreatinine ratio in spot urine sample was proposed by both the KDIGO and KDOQI guidelines as an alternative for follow-up of chronic kidney disease. However, in daily practice we have observed important differences in the values of 24-hour urine proteincreatinine ratio (24-UPC) and protein-creatinine ratio in spot urine sample (spot-UPC)

AIMS

- Analyze the concordance between 24-UPC and spot-UPC
- Analyze the correlation of both samples with 24-hour proteinuria

RESULTS

24-UPC has a tendency to be higher tan spot-UPC (Interquartile range of variants difference= 0,09-0,58)


- There is global correlation between both ratios - Wilcoxon p<0,001 - Spearman CI: 0,8876 PRO/CRO
- **Correlation between spot-UPC and 24-UPC**

METHODS

- Retrospective study
- Period March 2009-March 2014
- •N: 409 samples from 149 patients
- •Simultaneous 24-hour and spot urine samples at the same visit.


Exclusion criteria:

- •Patients with urine protein < 0,3 g/24h
- •Patients with monoclonal gammapathy
- •Renal transplant patients




Concordance between spot-UPC and 24-UPC

24-UPC cutoff	Area under the curve	Spot-UPC	Sensibility	Specificity
1	0,9184	0,8	84,87%	87,82%
3,5	0,9796	2,1	90,57%	90,17%



Correlation between 24-UPC and 24-hr proteinuria



Correlation between spot-UPC and 24-hr proteinuria

CONCLUSIONS

- There are differences between 24-hour urine-creatinine ratio and spot sample urine creatinine-ratio, with significant differences when >2,5 mg/mg
- 24-hour urine-creatinine ratio is more reliable for estimation of proteinuria in patients with proteinuric CKD.







